

MENU NAVIGATION

Lugitsch-Strasser GmbH

Program 121°C 00:20:00s
Phase STOPPED

Change Steril.temp/Steril.time

Servicemenü

Alarm overview

reset maintenace

Batches since last maintenace **600**
Total batches: **1200**
Maintenance: **600**

rest maintenace

select language

Deutsch

English

OK

calibrate sterilization temperature

calibrate drain temperature

calibrate pressure and external tank

reset maintenace

select language

Sterilizer number

Alarm overview

Technician
Password (xxxx)

Password
"0664"

Sterilizer number

Change sterilizer number and press OK

Yes

Sterilizer works automatical No

OK

Alarm overview

Alarm time	Event
18/05/16	20:14
18/05/16	20:22
18/05/16	20:34
18/05/16	20:34